

MONOZOOM MICROSCOPE JBM-MZ50

Digital Monozoom microscopes adopts apochromatic parallel optical imaging system and provide high resolution and sharp stereoscopic images, 4K HDMI digital camera with measurement function can work without PC. This microscope is designed for the applications of circuit board, electronics, semi-conductor and related industrial inspection and scientific research areas. Modularization product portfolio and excellent optical performance make it a better choice in these areas.

Magnification : 50x

Zoom objective : 0.7X-5X

Eyepiece : 0.4X

Zoom Ratio : 1:7

Working Distance: 100mm

Base Vertical post height : 240mm

Vertical post diameter : $\phi 32\text{mm}$

Microscope body holder : $\phi 50\text{mm}$

Focus travel : 50mm

Base size : 320x260x16mm

Material : Metal

(OPTIONAL) We Suggest the Camera According Application.



+91-9873571533
+91-9354741245



www.jbmicroscope.com
info@jbinstruments.com



J-9/36, S/F UNIT NO: 402-403,
RAJOURI GARDEN, DELHI, - 110027